

Introducing **v-gel**[®] for dogs



The first ever
supraglottic airway
devices for dogs
are **HERE** – and
they are even
better than you
imagined...



anaesthesia – safer with v-gel[®]



www.docsinnovent.com

What is a Supraglottic Airway Device (SGAD)?

Widely used in human medicine, a SGAD is an alternative to an endotracheal tube, used to supply anaesthetic gases and/or oxygen to anaesthetised patients. Designed to sit in the pharynx, they avoid many of the problems caused by ET tubes such as laryngeal/tracheal trauma and airway resistance. They are also quick and easy to place and are useful for emergency resuscitation. v-gel[®] are the first veterinary species-specific SGADs, currently used all over the world in cats and rabbits.



v-gel[®] advanced devices are suitable for a wide variety of procedures, including dental procedures



NEW to dog v-gel[®]

- **Gastric channel** – with insertion of oesophageal or gastric tube, allows for greater protection and management of reflux/regurgitation and aspiration
- **Single use** – human medicine standards of hygiene


SAME GREAT BENEFITS of v-gel[®] cat and v-gel[®] rabbit

- **Anatomical matching features** combined with a **soft material** to give a **high quality pressure seal**
- Perfectly suited to **ventilated** as well as spontaneously breathing patients
- **Avoids laryngeal and tracheal trauma**, which means **safer anaesthetics** and superb comfortable patient recoveries
- Fast, easy, safe and **stress-free** insertions
- **No post-operating coughing** or gagging
- **Low airway breathing resistance** due to the large airway channel within the device
- High quality pressure seal **restricting leakage of volatile anaesthetic agents**, thus improving health and safety for staff
- Suitable for a variety of patients including brachycephalics, and procedures including **dentals**

Key features of NEW dog v-gel[®] with integrated gastric channel

Gastric channel ingress – insert gastric tube to enable detection and treatment of reflux/regurgitation

Gastric channel egress – emerges beyond oesophageal seal, protecting pharynx from reflux/regurgitation





REFLUX/REGURGITATION IN DOGS DURING ANAESTHESIA

- Reflux and regurgitation during anaesthesia is much more common than many people think as they are not seeing it!
 - Reflux – 16–55% (Raptoulos and Galatos 1997)
 - Regurgitation – 6–15% (Wilson et al. 2006)
- Consequences of reflux regurgitation can be
 - Mild – oesophagitis in recovery period
 - Severe – oesophageal stricture, aspiration pneumonia
- Gastric channel in v-gel[®] can be used both to detect and treat (by flushing) reflux/regurgitation






v-gel[®] with gastric tube inserted



Sizing Guide

This table is intended to be a guide only. The right size v-gel[®] advanced will be determined by the size and shape of a dog's pharynx, which varies depending on the breed, size and head and neck shape of the dog.

Head shape – Brachycephalic	Lean weight of dog	Head shape – Mesocephalic	Lean weight of dog	Head shape – Dolichocephalic	Lean weight of dog	Size/ Product code
	6 – 10kg 13 – 22lb		4 – 7kg 9 – 15lb		2 – 4kg 4.5 – 9lb	D3 = D20003
	9 – 15kg 20 – 33lb		6 – 10kg 13 – 22lb		4 – 7kg 9 – 15lb	D4 = D20004
	14 – 20kg 30 – 45lb		9 – 15kg 20 – 33lb		6 – 10kg 13 – 22lb	D5 = D20005
	19 – 26kg 42 – 57lb		14 – 20kg 30 – 45lb		9 – 15kg 20 – 33lb	D6 = D20006

Also available



v-gel[®] advanced for cats and rabbits, proven technology
for safer anaesthesia

VetLube[®] Lubricant
a water-based sterile
lubricant spray for veterinary
use to assist in insertion or
placement of v-gel[®] advanced
supraglottic airway and other
medical devices.

Product code: D14002



**d-grip[®] circuit
support**
for support of
anaesthetic circuit
to prevent drag on
the v-gel[®] or other
airway devices.

**Product code:
D13001**

